



People.Health.Care.

VENOSAN[®] 5000


The magic of silver.



VENOSAN® 5000


Compression stockings with silver.

The premium product for those with higher expectations.




Pure silver-coated fibers

Presents the unique benefits of silver with its silver-coated X-Static® yarn directly against the skin.



Heat conductive

Equalizes temperature in both cold and warm climates by transporting moisture away from the skin.



Anti-static

Conducts static electricity charges away from everyday clothing with its silver-coated X-Static yarn.



Long-lasting

Maintains its silver properties up to 150 washing cycles.



Durable

Robust knit construction makes it a great option for active men and women.



Swiss comfort

Manufactured in Switzerland to the highest standards.



Available Styles




VENOSAN® 5000 medical grade compression stockings offer mid-containment of mild to moderate venous edema, early stages of lipedema and preventative support. Its leading benefit of X-Static silver-coated yarn provides the unique benefits of silver — making it a great option for those with sensitive skin. Additionally, it offers a soft touch and comfortable experience with its microfiber TACTEL® yarns that help regulate body temperature by rapidly wicking sweat and moisture away from the body.

Indications

- Chronic Venous Insufficiency (CVI) (Stage 1–2)
- Post vein surgery
- Mild skin conditions (e.g., irritation, pigmentation)
- Lipedema (Stage 1)

Options

Compression	20–30 mmHg, 30–40 mmHg, 40–50 mmHg
Toe Style	Open toe, Closed toe
Fit	Standard, Made to Measure
Colours	<div style="display: flex; gap: 10px;"> <div style="background-color: black; color: white; border-radius: 50%; padding: 5px 10px;">Black</div> <div style="background-color: #8B4513; color: white; border-radius: 50%; padding: 5px 10px;">Marokko</div> <div style="background-color: #A0522D; color: white; border-radius: 50%; padding: 5px 10px;">Mexico</div> <div style="background-color: #D2691E; color: white; border-radius: 50%; padding: 5px 10px;">Sand</div> </div>



A Silver threads visible on inner side of VENOSAN® 5000 under a microscope



People.Health.Care.


VENOSAN® 5000
Compression stockings with silver.

Measurements | Mesures de circonférence (cm)

	SMALL	MEDIUM	LARGE	X-LARGE
cB	20-23	23-26	26-29	29-31
cC	29-39	33-42	36-45	39-47
cD	28-38	32-41	35-44	38-46
cG	48-63	54-67	58-72	63-74

Length | Longueur (cm)

		short	long
AD/ADH	ℓD	35-40	37-42
AGH	ℓG	63-72	68-78
AGG/AT	ℓK	67-74	74-85



Material content | Composition
40% TACTEL® polyamide, 30% elastane, 20% cotton, 10% X-Static® — The Silver Fiber™



L&R CAN INC. · 2820 Argentia Road, Unit 2, Mississauga, ON L5N 8G4 · 888-836-6726 · Fax: 888-434-5347 · inquiries@ca.LRmed.com · www.Lohmann-Rauscher.ca

VENOSAN® is a registered trademark of SWISSLASTIC AG, ST. GALLEN, Switzerland — an L&R Company
TACTEL® fibre is a registered trademark of INVISTA and is used under license to NILIT for nylon fibre products.